INSTRUCTIONS FOR USE

AlproJet-D

Part 1 - Alkaline Cleaning | Step 1 of 2 | Aspiration systems

REF 3104 · 1 L bottle · REF 3101 · 5 L canister

ALPRO-AU-IFU-ALPROJET-D · Version 1.2 · Issue October 2026 · Next review October 2027



Part 1 alkaline cleaning for aspiration systems

Concentrated alkaline cleaner for routine cleaning and maintenance of dental aspiration and separation systems, including units with amalgam separators. AlproJet-D is Part 1 of the two-part cleaning protocol; AlproJet-W is Part 2 acidic cleaning.

Instructions for Use

1. **Prepare the solution -** Prepare 2 L total using warm water (40-50 C). Mix 50 mL AlproJet-D concentrate with 950 mL water twice, making two separate 1 L mixtures.
2. **Aspiration lines -** Use the first 1 L mixture and aspirate it through the suction lines.
3. **Spittoon -** Pour the second 1 L mixture through the spittoon.
4. **Contact time -** Allow 30-60 minutes for best cleaning results. If left in the system longer than recommended, AlproJet-D has no contact-time sensitivity when used as directed.
5. **Alternating schedule -** Use AlproJet-D and AlproJet-W on different days, following the schedule below.

Alternating-use examples

For centralised suction systems, with or without an amalgam separator: use AlproJet-D for one full week, followed by AlproJet-W for the next full week, alternating weekly.

For non-centralised suction systems: use AlproJet-D for 3 days and AlproJet-W for 2 days each week.

IMPORTANT: AlproJet-D and AlproJet-W must NEVER be used together or on the same day.

Delivery: REF 3104 - 1 L bottle · REF 3101 - 5 L canister.

Shelf life: 3 years.

System: Part 1 alkaline cleaning; Part 2 is AlproJet-W acidic cleaning.

Always follow the current Product Information and Safety Data Sheet and the equipment manufacturer's instructions. Biodegree supports compliance; it does not certify practice compliance.