INSTRUCTIONS FOR USE

Alpron

Daily Water Decontamination Agent

For dental units and other medical devices with external water bottle systems

REF 3184 · 1 L bottles · ARTG 226413

ALPRO-AU-IFU-ALPRON · Version 1.1 · Issue October 2026 · Next review October 2027



Daily decontamination for procedural water

Alpron is a liquid concentrate for daily decontamination and maintenance of procedural water in dental units and other medical devices supplied through external water bottles or canisters.

Daily use sequence

1. **Prepare a fresh 1% Alpron solution**
   For a 1 L bottle, mix 10 mL Alpron with 990 mL clean water. For larger containers, keep the same 1% ratio and scale the quantities proportionally.

Final volume	Alpron	Clean water
1 L	10 mL	990 mL
2 L	20 mL	1,980 mL
5 L	50 mL	4,950 mL

2. **Connect and use during the working day**
   Fill the external bottle with the prepared 1% solution, connect it to the dental unit, and use it as the procedural water supply throughout the day.

3. **Leave remaining solution overnight**
   If Alpron solution remains in the bottle at the end of the day, leave it in place overnight so it can continue acting as a decontaminant.

4. **Start each morning with fresh solution**
   The following morning, discharge the remaining overnight solution and prepare a new 1% Alpron solution. Repeat this daily throughout the working week.

5. **End-of-week bottle cleaning and sanitisation**
   At the end of the week, remove all remaining liquid and clean the bottle. Sanitise the bottle with BC-San 100 from ALPRO Medical according to its current Instructions for Use.

6. **On return, rinse and restart**
   Discharge the BC-San 100 solution, rinse the bottle thoroughly with clean water, then prepare and fill with a fresh 1% Alpron solution before reconnecting the bottle.

Integrated dosing systems

Where the dental unit has an integrated dosing device, fill the designated reservoir with undiluted Alpron. Have the dosing device adjusted by a qualified service engineer to deliver Alpron at the manufacturer's specified 1% dilution.

Always follow the current Alpron and BC-San 100 Instructions for Use and Safety Data Sheets. Application should be carried out by appropriately trained personnel.