PRODUCT SUMMARY

BIB forte

Instruments

REF 3024 · 1 L bottle · ARTG 167778 · Class IIb Medical Device

ALPRO-AU-SUM-BIBFORTE · Version 1.0 · Issue October 2026 · Next review October 2027



More Than Cleaning Before Sterilisation

Why only clean when you can **clean and disinfect in the same ultrasonic step?**

BIB forte combines effective pre-cleaning with **validated bactericidal disinfection**, helping reduce microbial contamination before instruments proceed through inspection, packaging and required sterilisation.

2% • 5 minutes in an ultrasonic device
Clean → Disinfect → Prepare for Sterilisation

Disinfection does not replace sterilisation where sterilisation is required. Always follow the reusable medical device manufacturer's IFU.

BIB forte — More Than Cleaning
Unlike an ordinary instrument cleaner, **BIB forte combines cleaning with validated medical-device disinfection**.

Key features

- **Clean + disinfect in one ultrasonic step.** At the validated **2% bactericidal concentration / 5-minute ultrasonic cycle**, BIB forte provides cleaning together with bactericidal disinfection.
- **Reduce microbial load before handling.** Instruments are not simply cleaner—the microbial load has been reduced before staff proceed with rinsing, inspection, drying and packaging.
- **Class IIb medical-device disinfectant.** BIB forte's Class IIb classification reflects its intended purpose as a **medical-device disinfectant**, rather than a cleaner alone.
- **Better preparation for sterilisation.**
  - **Ordinary cleaner:** Cleaning → Sterilisation
  - **BIB forte: Cleaning + Disinfection → Sterilisation**
- **Disinfection where appropriate.** For compatible devices where disinfection is the required endpoint, BIB forte may be used according to its validated claims and the device manufacturer's IFU.

Important: BIB forte does **not replace sterilisation** for critical or other reusable instruments that are required to be sterile.