INSTRUCTIONS FOR USE

Bilpron

Standby Waterline Protection

Protecting dental waterlines when water is not flowing

REF 3181 · ARTG 226413

ALPRO-AU-IFU-BILPRON · Version 1.1 · Issue October 2026 · Next review October 2027



Standby protection while your waterlines are at rest

Bilpron is a ready-to-use standby solution for weekends, holidays and other periods when water is not flowing. It remains inside the waterlines, coating internal surfaces and helping protect them from bacterial adhesion and re-colonisation until the unit is put back into service.

How to use Bilpron

1. **Prepare before closing**
   Before a weekend, holiday or other extended closure, fill the external water bottle with approximately 200-250 mL of undiluted Bilpron.
2. **Fill every waterline**
   Run Bilpron through all water outlets and waterlines until the solution is visibly coming through as blue water. This confirms that the lines are filled.
3. **Stop - do not rinse**
   Once blue Bilpron is visible at the outlets, stop the flow. Do not rinse. Leave Bilpron inside the waterlines so the internal surfaces remain covered during the standby period.
4. **Leave in place during closure**
   Bilpron remains effective in the waterlines during weekends and holidays. For extended shutdowns, it can remain in the system for up to 3 months, after which replenishment may be required.
5. **Flush on return**
   When the practice reopens, flush Bilpron out of all waterlines with clean water until the blue colour has completely disappeared and the water runs clear.
6. **Restart the daily Alpron protocol**
   Prepare the normal daily 1% Alpron solution, fill the water bottle and reconnect it to begin the new day's procedural-water protocol.

> **Important**
>
> Bilpron is a standby protection product, not a biofilm-removal treatment. It is intended to protect clean waterlines during periods without water flow. Use only as directed and follow the current Product Information, Safety Data Sheet and dental-unit manufacturer's instructions.