INSTRUCTIONS FOR USE

CleanWipes

Cleaning of Medical Device Surfaces

Class I Medical Device - Ready-to-Use Near-Neutral Detergent Wipes

REF 4153 · dispenser box with 70 wipes · REF 4154 · carton of 12 refill bags · ARTG 282859

ALPRO-AU-IFU-CLEANWIPES · Version 1.0 · Updated October 2026



Gentle on sensitive surfaces. Tough on organic contamination.

CleanWipes are ready-to-use, near-neutral detergent wipes for cleaning sensitive medical-device surfaces and dental equipment. They help dissolve and remove organic contamination, including blood, protein and tissue residues, while offering excellent material compatibility and minimal residue concerns.

The non-protein-fixing formulation helps lift contamination before it dries onto the surface, reduces microbial contamination and prepares the medical-device surface for effective subsequent disinfection where required.

CleanWipes are bacteriostatic and fungistatic. They are a cleaner, not a bactericidal or virucidal disinfectant.

Instructions for Use

1. **Ready to use**
   Remove one pre-moistened CleanWipe from the dispenser. No dilution or additional solution is required.
2. **Wipe the medical-device surface**
   Apply the wipe directly to the surface and distribute the cleaning solution evenly. Thoroughly moisten the surface, but avoid overwetting, pooling or excess liquid around openings and controls.
3. **Allow the cleaner to act on contamination**
   Where visible organic contamination is present, allow the moistened surface to remain wet briefly so blood, protein and tissue residues can soften and be lifted more effectively. Wipe until soiling is removed.
4. **Allow to dry naturally**
   After cleaning, allow the treated surface to dry naturally. Where disinfection is required, follow cleaning with an appropriate compatible medical-device disinfectant.
5. **Close the dispenser after use**
   Close the lid tightly after each wipe is removed to prevent drying. After opening, the wipes remain usable for at least 28 days when the lid is kept tightly closed.

Wipe characteristics

- **Large format**
  21 × 26 cm for greater surface coverage per wipe.
- **Viscose microfibre**
  Soft 'cotton touch' material for excellent solution absorption and even distribution across the treated surface.
- **Material compatibility**
  Near-neutral, alcohol-free cleaning chemistry for sensitive medical-device surfaces.
- **Environmental profile**
  Viscose wipe material offers a more environmentally compatible alternative to many generic polyester wipes.

> **Important**
>
> CleanWipes are intended for cleaning and decontamination of medical-device surfaces. They do not replace a required disinfection step. Always follow the current Product Information and Safety Data Sheet, and the medical-device manufacturer's cleaning and compatibility instructions.