INSTRUCTIONS FOR USE

MaxiWipes

Surfaces

REF 4612 · StarterSet: box + 250 dry wipes · REF 4613 · 4 × 250 refill · Regulated via the soaking solution

ALPRO-AU-IFU-MAXIWIPES · Version 1.0 · Issue October 2026 · Next review October 2027



Dry wipes for a compatible liquid.

Large, dry, lint-free non-woven wipes on a reel in the MaxiWipes-Box, for the cleaning and disinfection of the surfaces of medical devices in combination with PlastiSept eco or IC-100. The dry wipe has no efficacy of its own — the spectrum, contact time and Australian registration come from the solution it is soaked in.

Instructions for Use

1. **Qualified personnel only** — Usage only by qualified personnel. Always confirm the dry wipe is compatible with the chosen soaking solution.
2. **Open the box** — Remove the lid of the MaxiWipes Box. To refill, put the dry reel into the box.
3. **Soak the dry reel** — Soak the dry reel with the chosen disinfection solution: 1.5 L PlastiSept eco or 1.5 L IC-100 per box.
4. **Thread the first wipe and let it soak** — Pull the first wipe (from the centre of the reel) through the tear-off slot of the lid, close the box and allow the wipes to soak for 30 minutes. The wipes are then soaked sufficiently and ready for use.
5. **Label the box** — Complete the label of the MaxiWipes Box correctly: mark the checkbox for the soaking fluid used with a cross. Fixing the small label with the shelf life and lot number of the soaking fluid is recommended.
6. **Wipe — contact time of the solution** — Wipe the surface of the medical device. The spectrum of activity and the contact time are those stated in the Product Information of the soaking solution (PlastiSept eco, ARTG 164389 · IC-100, ARTG 238325).
7. **Keep the lid closed** — Always close the dispenser lid tightly. Once soaked, the wipes are usable for at least 28 days provided the lid is always closed tightly (manufacturer PI); discard the solution and wipes 1 month after preparation.

Field of application

1. Cleaning and disinfection of the surfaces of medical devices — the intended use follows the soaking solution.
2. Soaked with PlastiSept eco: especially suited to alcohol-sensitive surfaces of equipment such as treatment units (e.g. upholstery, head rests), suspension tables, the X-ray area (e.g. chin rests at OPG) and acrylic glass.
3. Soaked with IC-100: cleaning of medical-device surfaces, as stated in the IC-100 Product Information.

Wipe: tear-proof non-woven, “cotton touch” quality, 50 g/m², micro-pore technology, 17 × 26 cm. **Composition:** dry carrier with no active of its own. **Delivery:** REF 4612 — StarterSet: 1 dispenser box + 1 refill bag (250 dry wipes) · REF 4613 — carton of 4 refill bags (250 wipes each). **Safety:** the dry wipe carries no hazard of its own; once soaked, the Safety Data Sheet of the soaking disinfectant applies. **Disposal:** dispose of used wipes and residues per the soaking disinfectant’s SDS.

Always follow the current Product Information and Safety Data Sheet of the soaking solution, and the device manufacturer's cleaning and compatibility instructions. Biodegree supports compliance; it does not certify practice compliance.