ALPRO®
ALPRO MEDICAL GMBH

Biodegree Pty Ltd
www.biodegree.com.au
info@biodegree.com.au

HOW TO USE WL-CLEAN

EASY 3-SPRAY ROUTINE – CLEAN INSIDE IN UNDER 1 MINUTE

IMMEDIATE BENEFITS

LESS NOISE
Eliminates vibration & unwanted sounds

LESS HEAT
Reduces friction and internal stress

SMOOTHER RUN
Better performance, longer instrument life

1 SPRAY FIRST ~2 SEC



Insert adapter.
Spray ~ 2 sec.

2 WAIT 2 SEC



Wait 2 sec to allow solution to work.

3 SPRAY SECOND ~2 SEC



Spray again ~ 2 sec.

4 SPRAY THIRD ~2 SEC



Final spray ~ 2 sec.

5 DRY & ELIMINATE



Use WL-Blow (air gun) to blow out and dry channels.
No air gun?
Run lubrication cycle (2 times).

6 STERILISE



Sterilise in autoclave according to manufacturer's IFU.

WL-ADAPTERS – THE RIGHT SOLUTION FOR EVERY HANDPIECE

WL-Adapter 03/G* (Edelstahl/ stainless steel)	WL-Adapter 04/G* (Edelstahl/ stainless steel)	WL-Adapter 09/G* (Edelstahl/ stainless steel)	WL-Adapter 10/G* (Edelstahl/ stainless steel)	WL-Adapter 02/G* (Edelstahl/ stainless steel)
				
REF 3533	REF 3534	REF 3539	REF 3561	REF 3547
Sirona	KaVo	NSK	Morita	ISO

* Compliant with ISO 17664 & DIN 15804

WL-CLEAN 1 SPRAY



- Cleaning of spray channels & internal surfaces
- For all dental handpieces & instruments
- Medical device – Class I
- 500 ml aerosol can

NO CFC
FOR OUR ENVIRONMENT

MORE INFO: CONTACT

- ✓ Protect your instruments
- ✓ Reduce repairs & downtime
- ✓ Save money every day

info@biodegree.com.au