PRODUCT SUMMARY

# i3 TWO

i3 sterile filters · Shower

i3 TWO connect · i3 TWO direct · ordered through Biodegree · Not ARTG-registered

BD-AU-SUM-I3-TWO · Version 1.0 · Issue October 2026 · Next review October 2027



**For a compatible shower connection.**

The **i3 TWO** point-of-use sterile filter is used in medical settings to prevent infections caused by bacteria such as Legionella pneumophila, Pseudomonas aeruginosa and mycobacteria. Used as a shower head for the sterile filtration of piped water, it works where aerosols form and protects immunocompromised and immunosuppressed patients from waterborne pathogens.

Figures are manufacturer / EU data, not an Australian therapeutic claim; the i3 range is not ARTG-registered. Follow the manufacturer data sheet and Instructions for Use.

## Key features

- Bacteriostatic filter housing protects against retrograde contamination throughout the period of use.
- Stable flow rate throughout the 50-day period of use.
- Reliable sterile filtration with a high-performance flat membrane, nominal pore size 0.2 µm.
- Quick, tool-free installation — quick-release fastening (i3 TWO connect) or 1/2" screw thread (i3 TWO direct).
- CE-marked medical device to prevent infection; biocompatibility to ISO 10993-1.