PRODUCT SUMMARY

WL-clean

Instruments

REF 4150 · 4 × 500 g cans · ARTG 203023

ALPRO-AU-SUM-WLCLEAN · Version 1.0 · Issue October 2026 · Next review October 2027



Oiling Does Not Replace Internal Cleaning

Blood, saliva and fine debris can penetrate a handpiece's internal pathways, beyond the reach of external wiping.

WL-Clean provides ready-to-use, non-protein-fixing internal cleaning for compatible hollow instruments **before oiling and sterilisation**—addressing a critical step in reprocessing while helping maintain performance and **extend handpiece service life**.

Why Internal Cleaning Matters

Key features

- **Removes internal contamination** — helps eliminate blood, saliva, proteins and fine debris from internal pathways.
- **Prevents protein fixation and corrosion risk** — removes blood residues before they can become fixed to internal components and contribute to material deterioration.
- **Prepares for effective sterilisation** — supports appropriate pre-cleaning before subsequent maintenance and steam sterilisation.
- **Optimises lubrication** — cleaner internal pathways allow lubricating oil to reach and protect critical moving components more effectively.
- **Supports handpiece longevity** — helps preserve internal components, performance and overall service life.
- **Strengthens infection prevention and control** — adds an important internal-cleaning step to the handpiece reprocessing protocol, supporting safe preparation for subsequent patient use.